

Antibodies Product Datasheet

Product Name:	ZFP106 Antibody
Catalog Number:	R12-2415
Reactivity:	Human, Mouse, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZFP106. AA range:1833-1882
Concentration:	1 mg/ml
Gene Name:	ZFP106
Protein Name:	Zinc finger protein 106 homolog
Other Name:	ZFP106; SH3BP3; ZNF474; Zinc finger protein 106 homolog; Zfp-106; Zinc finger protein 474
Human Swiss Prot No.:	Q9H2Y7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ZFP106 Polyclonal Antibody detects endogenous levels of ZFP106 protein.
Storage and Stability:	-20°C/1 year