

Antibodies Product Datasheet

Product Name:	ZO-2 Antibody
Catalog Number:	R12-2417
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZO-2. AA range:1063-1112
Concentration:	1 mg/ml
Gene Name:	TJP2
Protein Name:	Tight junction protein ZO-2
Other Name:	TJP2; X104; ZO2; Tight junction protein ZO-2; Tight junction protein 2; Zona occludens protein 2; Zonula occludens protein 2
Human Swiss Prot No.:	Q9UDY2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ZO-2 Polyclonal Antibody detects endogenous levels of ZO-2 protein.
Storage and Stability:	-20°C/1 year