
Antibodies Product Datasheet

Product Name:	ACAD-10 Antibody
Catalog Number:	C14253
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ACAD10. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	ACAD10
Protein Name:	Acyl-CoA dehydrogenase family member 10
Other Name:	ACAD10; Acyl-CoA dehydrogenase family member 10; ACAD-10
Human Swiss Prot No.:	Q6JQN1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ACAD-10 Polyclonal Antibody detects endogenous levels of ACAD-10 protein.
Storage and Stability:	-20°C/1 year