
Antibodies Product Datasheet

Product Name:	A-FABP Monoclonal Antibody
Catalog Number:	N0013
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of A-FABP (aa61-121) expressed in E. Coli.
Concentration:	
Gene Name:	FABP4
Protein Name:	Fatty acid-binding protein, adipocyte
Other Name:	FABP4; Fatty acid-binding protein; adipocyte; Adipocyte lipid-binding protein; ALBP; Adipocyte-type fatty acid-binding protein; A-FABP; AFABP; Fatty acid-binding protein 4
Human Swiss Prot No.:	P15090
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	A-FABP Monoclonal Antibody detects endogenous levels of A-FABP protein.
Storage and Stability:	-20°C/1 year