

Antibodies Product Datasheet

Product Name:	ABCB5 Monoclonal Antibody(11A2)
Catalog Number:	N1339
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICC
Dilutions:	WB: 1:2000 IF 1:200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of ABCB5
Concentration:	
Gene Name:	ABCB5
Protein Name:	ATP-binding cassette sub-family B member 5
Other Name:	ABCB5; ATP-binding cassette sub-family B member 5; ABCB5 P-gp; P-glycoprotein ABCB5
Human Swiss Prot No.:	Q2M3G0
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous ABCB5 proteins.
Storage and Stability:	-20°C/1 year