

Antibodies Product Datasheet

Product Name:	14-3-3 eta Mouse mAb(5B2)
Catalog Number:	N2003
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Mouse,monoclonal
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human 14-3-3 eta
Concentration:	1 mg/ml
Gene Name:	YWHAH YWHA1
Protein Name:	14-3-3 eta
Other Name:	14-3-3 protein eta (Protein AS1)
Human Swiss Prot No.:	Q04917
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of 14-3-3 eta at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year