

Antibodies Product Datasheet

Product Name:	Myf-6 Antibody
Catalog Number:	C10457
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MYF6. AA range:116-165
Concentration:	1 mg/ml
Gene Name:	MYF6
Protein Name:	Myogenic factor 6
Other Name:	MYF6; BHLHC4; MRF4; Myogenic factor 6; Myf-6; Class C basic helix-loop-helix protein 4; bHLHc4; Muscle-specific regulatory factor 4
Human Swiss Prot No.:	P23409
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Myf-6 Polyclonal Antibody detects endogenous levels of Myf-6 protein.
Storage and Stability:	-20°C/1 year