

Antibodies Product Datasheet

Product Name:	PGP9.5 Rabbit Antibody(UCHL-1)
Catalog Number:	C42350
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein derived from human PGP9.5
Concentration:	1 mg/ml
Gene Name:	UCHL1
Protein Name:	PGP9.5
Other Name:	Ubiquitin carboxyl-terminal hydrolase isozyme L1 (UCH-L1) (EC 3.4.19.12) (EC 6.-.-) (Neuron cytoplasmic protein 9.5) (PGP 9.5) (PGP9.5) (Ubiquitin thioesterase L1); YM3801
Human Swiss Prot No.:	P09936
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human PGP9.5
Storage and Stability:	-20°C/1 year