

Antibodies Product Datasheet

Product Name:	CREB-1 Monoclonal Antibody
Catalog Number:	N0162
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF/ICELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human CREB-1 expressed in E. Coli.
Concentration:	
Gene Name:	CREB1
Protein Name:	Cyclic AMP-responsive element-binding protein 1
Other Name:	CREB1; Cyclic AMP-responsive element-binding protein 1; CREB-1; cAMP-responsive element-binding protein 1
Human Swiss Prot No.:	P16220
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	CREB-1 Monoclonal Antibody detects endogenous levels of CREB-1 protein.
Storage and Stability:	-20°C/1 year