

Antibodies Product Datasheet

Product Name:	RGAG1 Antibody
Catalog Number:	C11369
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RGAG1. AA range:1331-1380
Concentration:	1 mg/ml
Gene Name:	RGAG1
Protein Name:	Retrotransposon gag domain-containing protein 1
Other Name:	RGAG1; KIAA1318; Retrotransposon gag domain-containing protein 1; Tumor antigen BJ-HCC-23
Human Swiss Prot No.:	Q8NET4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RGAG1 Polyclonal Antibody detects endogenous levels of RGAG1 protein.
Storage and Stability:	-20°C/1 year