

## Antibodies Product Datasheet

|                               |  |
|-------------------------------|--|
| <b>Product Name:</b>          | EPHA7 Antibody   |
| <b>Catalog Number:</b>        | R12-2753   |
| <b>Reactivity:</b>            | Human, Mouse, Rat  |
| <b>Applications:</b>          | WB, ELISA  |
| <b>Dilutions:</b>             | EPHA7 Antibody   |
| <b>Modification:</b>          | Unmodified/Total   |
| <b>Source:</b>                | Rabbit Polyclonal  |
| <b>Purification:</b>          | The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.  |
| <b>Immunogen:</b>             | The antiserum was produced against synthesized peptide derived from human EPHA7.   |
| <b>Concentration:</b>         | 1 mg/ml  |
| <b>Gene Name:</b>             | EPHA7  |
| <b>Other Name:</b>            | Ephrin type-A receptor 7; Tyrosine-protein kinase receptor EHK-3; EPH homology kinase 3; Receptor protein-tyrosine kinase HEK11; EHK3; HEK11       |
| <b>Human Swiss Prot No.:</b>  | Q15375   |
| <b>Formulation:</b>           | Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| <b>Specificity:</b>           | EPHA7 Antibody antibody detects endogenous levels of EPHA7.  |
| <b>Storage and Stability:</b> | Stable at -20°C for at least 1 year.   |