

Antibodies Product Datasheet

Product Name:	GFPT2 Antibody
Catalog Number:	C50189
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GFPT2
Concentration:	1 mg/ml
Gene Name:	GFPT2
Protein Name:	Glucosamine--fructose-6-phosphate aminotransferase [isomerizing] 2 (EC 2.6.1.16) (D-fructose-6-phosphate amidotransferase 2) (Glutamine:fructose 6 phosphate amidotransferase 2) (GFAT 2) (GFAT2) (Hexosephosphate aminotransferase 2)
Other Name:	
Human Swiss Prot No.:	O94808
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GFPT2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year