

Antibodies Product Datasheet

Product Name:	DNAJC15 Antibody
Catalog Number:	C49875
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DNAJC15
Concentration:	1 mg/ml
Gene Name:	DNAJC15 DNAJD1 GIG22 HSD18
Protein Name:	DnaJ homolog subfamily C member 15 (Cell growth-inhibiting gene 22 protein) (Methylation-controlled J protein) (MCJ)
Other Name:	
Human Swiss Prot No.:	Q9Y5T4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DNAJC15 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year