

Antibodies Product Datasheet

Product Name:	Cytokeratin 16 (PT1870) mouse mAb
Catalog Number:	N1834
Reactivity:	Human, Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	WB 500-2000 IHC-p 1:100-500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse IgG
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Cytokeratin 16
Concentration:	
Gene Name:	KRT16 KRT16A
Protein Name:	Keratin, type I cytoskeletal 16 (Cytokeratin-16) (CK-16) (Keratin-16) (K16)
Other Name:	
Human Swiss Prot No.:	P08779
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Cytokeratin 16. Heat-induced epitope retrieval (HIER) Citrate buffer of pH6.0 was highly recommended as antigen repair method in paraffin section
Storage and Stability:	-20°C/1 year