

---

## Antibodies Product Datasheet

<b>Product Name:</b>	Src (Phospho-Tyr527) Antibody
<b>Catalog Number:</b>	A9474
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:1000-2000
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized phospho peptide around human Src (Tyr527)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	SRC SRC1
<b>Protein Name:</b>	Src (Tyr527)
<b>Other Name:</b>	Proto-oncogene tyrosine-protein kinase Src (EC 2.7.10.2) (Proto-oncogene c-Src) (pp60c-src) (p60-Src)
<b>Human Swiss Prot No.:</b>	P12931
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of Human Mouse Rat Src (phospho-Tyr527)
<b>Storage and Stability:</b>	-20°C/1 year