

## Antibodies Product Datasheet

<b>Product Name:</b>	Transferrin Monoclonal Antibody(5A4)
<b>Catalog Number:</b>	N1563
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB,IHC-p,IF(paraffin section),ELISA
<b>Dilutions:</b>	WB 1:2000-5000, IHC 1:100-200
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Protein
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	TF
<b>Protein Name:</b>	Serotransferrin (Transferrin) (Beta-1 metal-binding globulin) (Siderophilin)
<b>Other Name:</b>	TF; Serotransferrin; Transferrin; Beta-1 metal-binding globulin; Siderophilin
<b>Human Swiss Prot No.:</b>	P02787
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Transferrin protein detects endogenous levels of Transferrin
<b>Storage and Stability:</b>	-20°C/1 year