

Antibodies Product Datasheet

Product Name:	Cystatin C mouse Monoclonal Antibody(7F11)
Catalog Number:	N1608
Reactivity:	Human
Applications:	IF/ICC,WB,IHC-pELISA
Dilutions:	IF: 1:50-200 WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Cystatin C of CST3
Concentration:	1 mg/ml
Gene Name:	CST3
Protein Name:	CST3
Other Name:	CST3
Human Swiss Prot No.:	P01034
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cystatin C protein detects endogenous levels of CST3
Storage and Stability:	-20°C/1 year