

Antibodies Product Datasheet

Product Name:	GALNTL2 Antibody
Catalog Number:	C50405
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GALNTL2
Concentration:	1 mg/ml
Gene Name:	GALNTL2 UNQ770/PRO1564
Protein Name:	Polypeptide N-acetylgalactosaminyltransferase-like protein 2 (EC 2.4.1.41) (Polypeptide GalNAc transferase-like protein 2) (GalNAc-T-like protein 2) (pp-GaNTase-like protein 2) (Protein-UDP acetylgalactosaminyltransferase-like protein 2) (UDP-GalNAc:polyp
Other Name:	
Human Swiss Prot No.:	Q8N3T1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GALNTL2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year