

Antibodies Product Datasheet

Product Name:	TEF-5 Antibody
Catalog Number:	C20674
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TEAD3. AA range:211-260
Concentration:	1 mg/ml
Gene Name:	TEAD3
Protein Name:	Transcriptional enhancer factor TEF-5
Other Name:	TEAD3; TEAD5; TEF5; Transcriptional enhancer factor TEF-5; DTEF-1; TEA domain family member 3; TEAD-3
Human Swiss Prot No.:	Q99594
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TEF-5 Polyclonal Antibody detects endogenous levels of TEF-5 protein.
Storage and Stability:	-20°C/1 year