
Antibodies Product Datasheet

Product Name:	PI3 Kinase P85 Rabbit Antibody
Catalog Number:	C40137
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of PI3 Kinase P85 of PIK3R1
Concentration:	1 mg/ml
Gene Name:	PIK3R1
Protein Name:	PIK3R1
Other Name:	PIK3R1; YM3680
Human Swiss Prot No.:	P27986
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PI3 Kinase P85 protein detects endogenous levels of PIK3R1
Storage and Stability:	-20°C/1 year