

Antibodies Product Datasheet

Product Name:	Aurora A Antibody
Catalog Number:	C42088
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Aurora A
Concentration:	1 mg/ml
Gene Name:	AURKA AIK AIRK1 ARK1 AURA AYK1 BTAK IAK1 STK15 STK6
Protein Name:	Aurora A
Other Name:	Aurora kinase A (EC 2.7.11.1) (Aurora 2) (Aurora/IPL1-related kinase 1) (ARK-1) (Aurora-related kinase 1) (hARK1) (Breast tumor-amplified kinase) (Serine/threonine-protein kinase 15) (Serine/threonine-protein kinase 6) (Serine/threonine-protein kinase aur
Human Swiss Prot No.:	O14965
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.20% sodium azide.
Specificity:	This antibody detects endogenous levels of Aurora A at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year