

Antibodies Product Datasheet

Product Name:	ITGA3 (heavy chain, Cleaved-Phe33) Antibody
Catalog Number:	L0391
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ITGA3 (heavy chain, Cleaved-Phe33)
Concentration:	1 mg/ml
Gene Name:	ITGA3 MSK18
Protein Name:	ITGA3 (heavy chain, Cleaved-Phe33)
Other Name:	Integrin alpha-3 (CD49 antigen-like family member C) (FRP-2) (Galactoprotein B3) (GAPB3) (VLA-3 subunit alpha) (CD antigen CD49c) [Cleaved into: Integrin alpha-3 heavy chain; Integrin alpha-3 light chain]
Human Swiss Prot No.:	P26006
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat ITGA3 (heavy chain, Cleaved-Phe33, protein was cleaved amino acid sequence between 32-33)
Storage and Stability:	-20°C/1 year