

Antibodies Product Datasheet

Product Name:	Cbl Monoclonal Antibody
Catalog Number:	N0095
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICCFCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Cbl expressed in E. Coli.
Concentration:	
Gene Name:	CBL
Protein Name:	E3 ubiquitin-protein ligase CBL
Other Name:	CBL; CBL2; RNF55; E3 ubiquitin-protein ligase CBL; Casitas B-lineage lymphoma proto-oncogene; Proto-oncogene c-Cbl; RING finger protein 55; Signal transduction protein CBL
Human Swiss Prot No.:	P22681
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	Cbl Monoclonal Antibody detects endogenous levels of Cbl protein.
Storage and Stability:	-20°C/1 year