

Antibodies Product Datasheet

Product Name:	CYP17A1 Antibody
Catalog Number:	C12247
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cytochrome P450 17A1. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	CYP17A1
Protein Name:	Steroid 17-alpha-hydroxylase/17,20 lyase
Other Name:	CYP17A1; CYP17; S17AH; Steroid 17-alpha-hydroxylase/17; 20 lyase; CYPXVII; Cytochrome P450 17A1; Cytochrome P450-C17; Cytochrome P450c17; Steroid 17-alpha-monooxygenase
Human Swiss Prot No.:	P05093
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CYP17A1 Polyclonal Antibody detects endogenous levels of CYP17A1 protein.
Storage and Stability:	-20°C/1 year