

Antibodies Product Datasheet

Product Name:	STAM2 (Phospho-Tyr192) Antibody
Catalog Number:	A1237
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human STAM2 around the phosphorylation site of Tyr192. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	STAM2
Protein Name:	Signal transducing adapter molecule 2
Other Name:	STAM2; HBP; Signal transducing adapter molecule 2; STAM-2; Hrs-binding protein
Human Swiss Prot No.:	O75886
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-STAM2 (Y192) Polyclonal Antibody detects endogenous levels of STAM2 protein only when phosphorylated at Y192.
Storage and Stability:	-20°C/1 year