

Antibodies Product Datasheet

Product Name:	CD88 Antibody
Catalog Number:	B0853
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CD88/C5aR. AA range:301-350
Concentration:	1 mg/ml
Gene Name:	C5AR1
Protein Name:	C5a anaphylatoxin chemotactic receptor
Other Name:	C5AR1; C5AR; C5R1; C5a anaphylatoxin chemotactic receptor; C5a-R; C5aR; CD antigen CD88
Human Swiss Prot No.:	P21730
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD88 Polyclonal Antibody detects endogenous levels of CD88 protein.
Storage and Stability:	-20°C/1 year