

## Antibodies Product Datasheet

<b>Product Name:</b>	Ub (Acetyl-Lys11) Antibody
<b>Catalog Number:</b>	D12234
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB 1:1000-2000 ELISA 1:5000-20000
<b>Modification:</b>	Acetyl Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human Ub (Acetyl Lys11)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	UBA52 UBCEP2
<b>Protein Name:</b>	Ub (Acetyl Lys11)
<b>Other Name:</b>	Ubiquitin-60S ribosomal protein L40 (CEP52) (Ubiquitin A-52 residue ribosomal protein fusion product 1) [Cleaved into: Ubiquitin; 60S ribosomal protein L40]
<b>Human Swiss Prot No.:</b>	P62987/P62979/P0CG47/P0CG48
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of Human, Mouse, Rat Ub (Acetyl Lys11)
<b>Storage and Stability:</b>	-20°C/1 year