

Antibodies Product Datasheet

Product Name:	Cyclophilin E Antibody
Catalog Number:	C50912
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Cyclophilin E
Concentration:	1 mg/ml
Gene Name:	PPIE CYP33
Protein Name:	Peptidyl-prolyl cis-trans isomerase E (PPIase E) (EC 5.2.1.8) (Cyclophilin E) (Cyclophilin-33) (Rotamase E)
Other Name:	
Human Swiss Prot No.:	Q9UNP9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Cyclophilin E at Human, Mouse
Storage and Stability:	Store at -20° C/1 year