

Antibodies Product Datasheet

Product Name:	WAVE1 Antibody
Catalog Number:	B0598
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human WAVE1. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	WASF1
Protein Name:	Wiskott-Aldrich syndrome protein family member 1
Other Name: Human Swiss Prot No.:	WASF1; KIAA0269; SCAR1; WAVE1; Wiskott-Aldrich syndrome protein family member 1; WASP family protein member 1; Protein WAVE-1; Verprolin homology domain-containing protein 1 Q92558
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	WAVE1 Polyclonal Antibody detects endogenous levels of WAVE1 protein.
Storage and Stability:	-20°C/1 year