

Antibodies Product Datasheet

Product Name:	Spectrin β II Antibody
Catalog Number:	C18903
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SPTBN1. AA range:651-700
Concentration:	1 mg/ml
Gene Name:	SPTBN1
Protein Name:	Spectrin beta chain non-erythrocytic 1
Other Name:	SPTBN1; SPTB2; Spectrin beta chain; non-erythrocytic 1; Beta-II spectrin; Fodrin beta chain; Spectrin, non-erythroid beta chain 1
Human Swiss Prot No.:	Q01082
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Spectrin β II Polyclonal Antibody detects endogenous levels of Spectrin β II protein.
Storage and Stability:	-20°C/1 year