

Antibodies Product Datasheet

Product Name:	POU3F4 Antibody
Catalog Number:	C50049
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human POU3F4
Concentration:	1 mg/ml
Gene Name:	POU3F4 BRN4 OTF9
Protein Name:	POU domain, class 3, transcription factor 4 (Brain-specific homeobox/POU domain protein 4) (Brain-4) (Brn-4) (Octamer-binding protein 9) (Oct-9) (Octamer-binding transcription factor 9) (OTF-9)
Other Name:	
Human Swiss Prot No.:	P49335
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of POU3F4 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year