

Antibodies Product Datasheet

Product Name:	INTU Antibody
Catalog Number:	C50606
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human INTU
Concentration:	1 mg/ml
Gene Name:	INTU KIAA1284 PDZD6 PDZK6
Protein Name:	Protein inturned (Inturned planar cell polarity effector homolog) (PDZ domain-containing protein 6)
Other Name:	
Human Swiss Prot No.:	Q9ULD6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of INTU at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year