
Antibodies Product Datasheet

Product Name:	CD37 Monoclonal Antibody
Catalog Number:	N0120
Reactivity:	Human, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of CD37 expressed in E. Coli.
Concentration:	
Gene Name:	CD37
Protein Name:	Leukocyte antigen CD37
Other Name:	CD37; TSPAN26; Leukocyte antigen CD37; Tetraspanin-26; Tspan-26; CD antigen CD37
Human Swiss Prot No.:	P11049
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	CD37 Monoclonal Antibody detects endogenous levels of CD37 protein.
Storage and Stability:	-20°C/1 year