

---

## Antibodies Product Datasheet

<b>Product Name:</b>	Bc10 Antibody
<b>Catalog Number:</b>	C10782
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IF/ICC, ELISA
<b>Dilutions:</b>	WB 1:500-2000 Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human BLCAP. AA range:38-87
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	BLCAP
<b>Protein Name:</b>	Bladder cancer-associated protein
<b>Other Name:</b>	BLCAP; BC10; Bladder cancer-associated protein; Bladder cancer 10 kDa protein; Bc10
<b>Human Swiss Prot No.:</b>	P62952
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Bc10 Polyclonal Antibody detects endogenous levels of Bc10 protein.
<b>Storage and Stability:</b>	-20°C/1 year