

Antibodies Product Datasheet

Product Name:	CTDSP2 Antibody
Catalog Number:	C50099
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CTDSP2
Concentration:	1 mg/ml
Gene Name:	CTDSP2 NIF2 OS4 SCP2
Protein Name:	Carboxy-terminal domain RNA polymerase II polypeptide A small phosphatase 2 (EC 3.1.3.16) (Nuclear LIM interactor-interacting factor 2) (NLI-interacting factor 2) (Protein OS-4) (Small C-terminal domain phosphatase 2) (Small CTD phosphatase 2) (SCP2)
Other Name:	
Human Swiss Prot No.:	O14595
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CTDSP2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year