

Antibodies Product Datasheet

Product Name:	ACTR-IB Antibody
Catalog Number:	C10576
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ACV1B. AA range:73-122
Concentration:	1 mg/ml
Gene Name:	ACVR1B
Protein Name:	Activin receptor type-1B
Other Name:	ACVR1B; ACVRLK4; ALK4; Activin receptor type-1B; Activin receptor type IB; ACTR-IB; Activin receptor-like kinase 4; ALK-4; Serine/threonine-protein kinase receptor R2; SKR2
Human Swiss Prot No.:	P36896
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ACTR-IB Polyclonal Antibody detects endogenous levels of ACTR-IB protein.
Storage and Stability:	-20°C/1 year