
Antibodies Product Datasheet

Product Name:	Cdc6 Antibody
Catalog Number:	B0413
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDC6. AA range:20-69
Concentration:	1 mg/ml
Gene Name:	CDC6
Protein Name:	Cell division control protein 6 homolog
Other Name:	CDC6; CDC18L; Cell division control protein 6 homolog; CDC6-related protein; Cdc18-related protein; HsCdc18; p62(cdc6); HsCDC6
Human Swiss Prot No.:	Q99741
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cdc6 Polyclonal Antibody detects endogenous levels of Cdc6 protein.
Storage and Stability:	-20°C/1 year