

Antibodies Product Datasheet

Product Name:	Sec16A Antibody
Catalog Number:	C18482
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SEC16A. AA range:1013-1062
Concentration:	1 mg/ml
Gene Name:	SEC16A
Protein Name:	Protein transport protein Sec16A
Other Name:	SEC16A; KIAA0310; SEC16; SEC16L; Protein transport protein Sec16A; SEC16 homolog A
Human Swiss Prot No.:	O15027
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Sec16A Polyclonal Antibody detects endogenous levels of Sec16A protein.
Storage and Stability:	-20°C/1 year