
Antibodies Product Datasheet

Product Name:	p16 Antibody
Catalog Number:	C0285
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human p16 INK. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	CDKN2A
Protein Name:	Cyclin-dependent kinase inhibitor 2A isoforms 1/2/3
Other Name:	CDKN2A; CDKN2; MTS1; Cyclin-dependent kinase inhibitor 2A; isoforms 1/2/3; Cyclin-dependent kinase 4 inhibitor A; CDK4I; Multiple tumor suppressor 1; MTS-1; p16-INK4a; p16-INK4; p16INK4A
Human Swiss Prot No.:	P42771
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	p16 Polyclonal Antibody detects endogenous levels of p16 protein.
Storage and Stability:	-20°C/1 year