

Antibodies Product Datasheet

Product Name:	α -tubulin Monoclonal Antibody(8F11), HRP Conjugated
Catalog Number:	N1241
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF/ICCIP
Dilutions:	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB: 1:5000-10000, IHC 1:50-300.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse IgG
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	
Concentration:	1mg/ml
Gene Name:	TUBA1A
Protein Name:	Tubulin alpha-1A chain
Other Name:	TUBA1A
Human Swiss Prot No.:	Q71U36/P68363
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Specificity:	α -tubulin Monoclonal Antibody(8F11) HRP Conjugated, specially designed for your Western blot analysis.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.