

## Antibodies Product Datasheet

<b>Product Name:</b>	IL-17RB Antibody
<b>Catalog Number:</b>	C30894
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB 1:500-2000, ELISA 1:10000-20000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from IL-17RB at AA range: 251-300
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	IL17RB
<b>Protein Name:</b>	IL-17RB
<b>Other Name:</b>	Interleukin-17 receptor B (IL-17 receptor B) (IL-17RB) (Cytokine receptor-like 4) (IL-17 receptor homolog 1) (IL-17Rh1) (IL17Rh1) (Interleukin-17B receptor) (IL-17B receptor)
<b>Human Swiss Prot No.:</b>	Q9NRM6
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	IL-17RB Polyclonal Antibody detects endogenous levels of IL-17RB
<b>Storage and Stability:</b>	-20°C/1 year