
Antibodies Product Datasheet

Product Name:	MYL7 Antibody
Catalog Number:	C20920
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MYL7. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	MYL7
Protein Name:	Myosin regulatory light chain 2 atrial isoform
Other Name:	MYL7; MYL2A; MYLC2A; Myosin regulatory light chain 2; atrial isoform; MLC-2a; MLC2a; Myosin light chain 2a; Myosin regulatory light chain 7
Human Swiss Prot No.:	Q01449
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MYL7 Polyclonal Antibody detects endogenous levels of MYL7 protein.
Storage and Stability:	-20°C/1 year