

Antibodies Product Datasheet

Product Name:	Anti-Human BOB.1 (ABT-BOB1) Mouse Monoclonal Antibody
Catalog Number:	V0009
Clone ID:	ABT-BOB1
Reactivity:	Human, Mouse, Rat
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	POU2AF1
Gene Name:	Q16633
Other Name:	B cell Oct binding protein 1, B cell specific coactivator OBF 1, B cell specific coactivator OBF1, B-cell-specific coactivator OBF-1, BOB 1, BOB-1, OBF 1, OBF1, OBF1_HUMAN, OCA B, OCA-B, OCAB, OCT binding factor 1, OCT-binding factor 1, POU class 2 associating factor 1, POU domain class 2 associating factor 1, POU domain class 2-associating factor 1, Pou2af1
Isotype	IgG2a, Kappa
Localization:	Tonsil
Positive Tissue Control:	Nuclear, Cytoplasmic
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human BOB.1 protein.
Storage and Stability:	-20°C/1 year