47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145

www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: ICAM-1 Antibody

Catalog Number: B7113

Reactivity: Human, Rat, Mouse

Applications: WB, IHC-p, IF(paraffin section), ELISA

Dilutions: Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other

applications.

Modification: Unmodified/Total

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human ICAM-1. AA range:479-528

Concentration: 1 mg/ml

Gene Name: ICAM1

Protein Name: Intercellular adhesion molecule 1

Other Name: ICAM1; Intercellular adhesion molecule 1; ICAM-1; Major group

rhinovirus receptor; CD antigen CD54

Human Swiss Prot No.: P05362

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: ICAM-1 Polyclonal Antibody detects endogenous levels of ICAM-1

protein.

Storage and Stability: -20°C/1 year