

Antibodies Product Datasheet

Product Name:	VPS37B Antibody
Catalog Number:	C50664
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human VPS37B
Concentration:	1 mg/ml
Gene Name:	VPS37B
Protein Name:	Vacuolar protein sorting-associated protein 37B (hVps37B) (ESCRT-I complex subunit VPS37B)
Other Name:	
Human Swiss Prot No.:	Q9H9H4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of VPS37B at Human, Mouse
Storage and Stability:	Store at -20° C/1 year