

Antibodies Product Datasheet

Product Name:	ATP5A Antibody
Catalog Number:	C14602
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATP5A1. AA range:201-250
Concentration:	1 mg/ml
Gene Name:	ATP5A1
Protein Name:	ATP synthase subunit alpha mitochondrial
Other Name:	ATP5A1; ATP5A; ATP5AL2; ATPM; ATP synthase subunit alpha; mitochondrial
Human Swiss Prot No.:	P25705
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ATP5A Polyclonal Antibody detects endogenous levels of ATP5A protein.
Storage and Stability:	-20°C/1 year