

Antibodies Product Datasheet

Product Name:	PBF Antibody
Catalog Number:	C10723
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PTTG. AA range:116-165
Concentration:	1 mg/ml
Gene Name:	PTTG1IP
Protein Name:	Pituitary tumor-transforming gene 1 protein-interacting protein
Other Name:	PTTG1IP; C21orf1; C21orf3; Pituitary tumor-transforming gene 1 protein-interacting protein; Pituitary tumor-transforming gene protein-binding factor; PBF; PTTG-binding factor
Human Swiss Prot No.:	P53801
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PBF Polyclonal Antibody detects endogenous levels of PBF protein.
Storage and Stability:	-20°C/1 year