

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | DR5 mouse mAb |
| Catalog Number: | N0924 |
| Reactivity: | Human |
| Applications: | WB |
| Dilutions: | wb 1:2000 |
| Modification: | Unmodified/ Total |
| Source: | Monoclonal Mouse |
| Purification: | The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | Purified recombinant human DR5 protein fragments expressed in E.coli. |
| Concentration: | 1 mg/ml |
| Gene Name: | dr5 |
| Protein Name: | |
| Other Name: | Fas like protein ;Apoptosis inducing protein TRICK2A/2B ;Apoptosis inducing receptor TRAIL R2 ;CD 262 ;CD262 ;CD262 antigen ;Cytotoxic TRAIL receptor 2 ;Death domain containing receptor for TRAIL/Apo 2L ;Death domain containing receptor for TRAIL/Apo2L ;D |
| Human Swiss Prot No.: | O14763 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | This antibody detects endogenous levels of DR5 and does not cross-react with related proteins. |
| Storage and Stability: | -20°C/1 year |