

Antibodies Product Datasheet

Product Name:	ENaC γ Antibody
Catalog Number:	R12-2131
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ENaC gamma. AA range:132-181
Concentration:	1 mg/ml
Gene Name:	SCNN1G
Protein Name:	Amiloride-sensitive sodium channel subunit gamma
Other Name:	SCNN1G; Amiloride-sensitive sodium channel subunit gamma; Epithelial Na(+) channel subunit gamma; ENaCG; Gamma-ENaC; Gamma-NaCH; Nonvoltage-gated sodium channel 1 subunit gamma; SCNEG
Human Swiss Prot No.:	P51170
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ENaC γ Polyclonal Antibody detects endogenous levels of ENaC γ protein.
Storage and Stability:	-20°C/1 year