

Antibodies Product Datasheet

Product Name:	Casein Kinase II β (Phospho-Ser209) Antibody
Catalog Number:	A0869
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CKII-beta around the phosphorylation site of Ser209. AA range:166-215
Concentration:	1 mg/ml
Gene Name:	CSNK2B
Protein Name:	Casein kinase II subunit beta
Other Name:	CSNK2B; CK2N; G5A; Casein kinase II subunit beta; CK II beta; Phosvitin; Protein G5a
Human Swiss Prot No.:	P67870
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Casein Kinase II β (S209) Polyclonal Antibody detects endogenous levels of Casein Kinase II β protein only when phosphorylated at S209.
Storage and Stability:	-20°C/1 year